

MIC-Strip Colistin

EUCAST recommendation

Susceptibility testing of bacteria by the broth microdilution method to determine the minimum inhibitory concentration (MIC)

- ▶ Several methods under the spotlight and evaluated by EUCAST
- ▶ Microdilution method was ranked first
- ▶ Reference method for susceptibility testing of Enterobacteriaceae, *Pseudomonas aeruginosa* and *Acinetobacter* spp. related to Colistin according to EN ISO 20776-1
- ▶ Easy visual reading



According to the EUCAST and CLSI recommendations of March 2016 and the study of November 2016

1	2	3	4	5	6	7	8	9	10	11	12
GC	COL	COL	COL	COL	COL	COL	COL	COL	COL	COL	COL
	0,0625	0,125	0,25	0,5	1	2	4	8	16	32	64